



### MAIN APPLICATIONS

- Emergency Dept.
- General Surgery
- Gynecology
- Liposuction
- Neurosurgery
- Obstetrics
- Operating Theatre
- Dental Practice

NEW HOSPIVAC 350 has been designed for professional aspiration of bodily fluids, tissues or bones of the patient during or after surgery. The state-of-the-art oilless pump guarantees high performances with excellent suction capacities and max vacuum reachable within a few seconds, with no need for maintenance. A clear dashboard along with a full range of accessories makes it the ideal device for surgical suction and liposuction.

AVAILABLE MODELS AND ACCESSORIES INCLUDED	BASIC	FS	FULL
Liquid Collection Jar with overflow valve system (different options, see below)	o	o	o
8x14mm Silicone Tubes, length 150cm (autoclavable) + Conical Connector ø 10-11-12mm	o	o	o
Antibacterial & Hydrophobic Filter (single-patient)	o	o	o
Power Cord with Schuko plug	o	o	o
Footswitch (may be ordered at a later time)		o	o
Electronic Change-Over System from jar to jar			o

### TECHNICAL FEATURES

Motor	Oiless and maintenance-free piston pump
Power Feeding	220/230V-50/60Hz - other configurations on request
ISO 10079-1 Classification	HIGH VACUUM / HIGH FLOW
Max Vacuum (adjustable)	-0.90 bar -90 kPa -675 mmHg
Max free air flow rate	60 l/min
Noise Level	51,7 dBA
Duty cycle	Non-stop operation
Weight	kg 13
Size	cm 46 x 85 x 42
Years of Warranty	2
Shipping carton	1
Place of Manufacturing	Italy

Italian  
Medical  
Touch  
MADE IN ITALY

### NEW HOSPIVAC CONFIGURATIONS AND JARS INTERCHANGE

HOSPIVAC series offer a great versatility to users when it comes to change type of jars for liquid collection, regardless in which configuration the suction unit has been supplied in origin. A full line of accessories complete the range. The tables in the next pages provide all the information for ordering.

CE 0123