



CPAP GGM NP SERIES

INFANT BUBBLE CPAP SYSTEM

Use in Neonate CPAP and Bubble CPAP System

GGM NP series Infant Bubble CPAP System, complete system including servo control humidifier, full ranges of nasal prong, bonnet, CPAP generator, heater wire circuit, pressure release safety valve.

Complete system for application of continuous positive airway pressure to neonates and infants by nasal prongs (NCPAP) administered with a heater wire circuit used in conjunction with a continuous flow source infant ventilator or a suitably equipped multipurpose ventilator.



Cart, Trolley

Air/O2 Blender

INFANT BUBBLE CPAP SYSTEM / 嬰兒用氣泡式正壓通氣組

REF / 型號	Description / 內容	Pack / 包裝
VN-BC01 Bubble CPAP SYSTEM	1 x Bubble CPAP circuit w/chamber 1 x Safety valve kit 1 x CPAP generator	20 Sets/ctn
VN-BC02 Bubble CPAP SYSTEM W/STARTING KIT	1 x VH-2600A Servo control humidifier w/heater wire adaptor & temp probe 3 x Bubble CPAP circuit w/chamber 3 x Safety valve kit 3 x CPAP Generator 3 x Nasal prong kit (All size:0~5) 3 x Bonnet (All size: S, M, L)	3 Sets/ctn

Air/O2 Blender / 空氣氧氣混合器

REF / 型號	Description / 內容	Pack / 包裝
VN-BD01	Air/O2 Blender 2-18 LPM 1 x Air pressure line 1 x O2 pressure line 1 x C-Clamp	8 Sets/ctn
VN-BD02	Air/O2 Blender 5-50 LPM 1 x Air pressure line 1 x O2 pressure line 1 x C-Clamp	8 Sets/ctn

Cart, Trolley / 推車

REF / 型號	Description / 內容	Pack / 包裝
VN-BD04	With support arm & water bag arm	1Set/ctn

Nasal Prong Kit / 嬰兒鼻管組

REF / 型號	Description / 內容	Pack / 包裝
VN-0001	#0 (7.5fr) < 700g	10 Sets/box
VN-0002	#1 (10fr) < 700~1250g	10 Sets/box
VN-0003	#2 (12.5fr) < 1250~2000g	10 Sets/box
VN-0004	#3 (15fr) < 2000~3000g	10 Sets/box
VN-0005	#4 (16fr) < 3000g	10 Sets/box
VN-0006	#5 (20fr) < 1~2 years	10 Sets/box

Head Cap, Bonnet / 頭套

REF / 型號	Description / 內容	Pack / 包裝
VN-HC01	#S Size 11x13cm	10 Sets/box
VN-HC02	#M Size 12x14.5cm	10 Sets/box
VN-HC03	#L Size 12x15cm	10 Sets/box

Head strap / 頭帶

REF / 型號	Head circumference / 頭圍尺寸	Pack / 包裝
VN-HC11	#S Size 26cm~20cm	10 Sets/box
VN-HC12	#M Size 34cm~26cm	10 Sets/box
VN-HC13	#L Size 44cm~34cm	10 Sets/box



BUBBLE CPAP

Order Information



CART, TROLLEY

REF: NO>:VN-BD04

Description: With support arm & water bag arm

Pack: 1 Set/ctn

The whole system completes with a servo control humidifier, full ranges of nasal prong, CPAP generator, heater wire circuits and pressure release safety valve.



INFANT BUBBLE CPAP SYSTEM

REF	Description	Pack
VN-BC01	1 x Bubble CPAP circuit w/chamber 1 x Safety valve kit 1 x CPAP generator	20 Sets/ctn
VN-BC02	1 x VH-2600A Servo control humidifier w/heater wire adaptor & temp probe 3 x Bubble CPAP circuit w/chamber 3 x Safety valve kit 3 x CPAP Generator 3 x Nasal prong kit (All size: 0~5) 3 x Bonnet (All size: S, M, L)	5 Sets/ctn
VN-BC03	1 x Bubble CPAP circuit (Spiral)/chamber 1 x Safety valve kit 1 x CPAP generator	20 Sets/ctn
VN-BC04	1 x Bubble CPAP circuit (w/o heater wire) 1 x CO2 gas sampling line 1 x Bonnet 1 x Nasal prong kit	20 Sets/ctn

NASAL PRONG KIT

REF	Description	Pack
VN-0000	#00 (5fr) < 500g	10 Sets/box
VN-0001	#0 (7.5fr) < 700g	10 Sets/box
VN-0002	#1 (10fr) < 700~1250g	10 Sets/box
VN-0003	#2 (12.5fr) < 1250~2000g	10 Sets/box
VN-0004	#3 (15fr) < 2000~3000g	10 Sets/box
VN-0005	#4 (16fr) < 3000g	10 Sets/box
VN-0006	#5 (20fr) < 1~2 (years old)	10 Sets/box

HEAD CAP, BONNET

REF	Description	Pack
VN-HC01	#S Size 11x13cm	10 Sets/box
VN-HC02	#M Size 12x14.5cm	10 Sets/box
VN-HC03	#L Size 12x15cm	10 Sets/box

HEAD STRAP

REF	Head circumference	Pack
VN-HC11	#S Size 26cm~20cm	10 Sets/box
VN-HC12	#M Size 34cm~26cm	10 Sets/box
VN-HC13	#L Size 44cm~34cm	10 Sets/box



AIR/O2 BLENDER

REF	Description	Pack
VN-BD01	2-18 LPM 1 x Air pressure line 1 x O2 pressure line 1 x C-Clamp	8 Sets/ctn
VN-BD02	5-50 LPM 1 x Air pressure line 1 x O2 pressure line 1 x C-Clamp	8 Sets/ctn



GREAT GROUP MEDICAL CO., LTD.

No.168, Xingong 2nd Rd., Tianzhong Township,
Changhua County, Taiwan

Tel:886-4-875-8181 Fax:886-4-875-6161



info@greatgroup.com.tw

www.greatgroup.com.tw