

**Appearance:**

Stream-line design

Water drainage design

Compatible with other standard chambers

流線、排水設計，並相容於其他加濕瓶

**GENERAL PURPOSE HUMIDIFIER / 泛用型加熱器**

REF / 型號	VH-2100-115 / VH-2100-230
Dimensions / 規格	W135 x D170 x H156 mm
Weight / 重量	1.5 kg (w/o Chamber fitted)
Electrical Power / 電源	(1) Heater plate capacity : 150W max (2) Temperature Power cutout at $93^{\circ}\text{C} \pm 6^{\circ}\text{C}$ (3) $230 \pm 23\text{V}$ , 0.7A max. at 230V (4) $115 \pm 12\text{V}$ , 1.3 A max. at 115V
Heater Plate Overheat / 過熱保護	Cutoff: $93 \pm 5^{\circ}\text{C}$
Temperature Setting / 溫度設定	(1) From step 1 to step 9 (2) Range $45 \sim 77^{\circ}\text{C}$
Standards / 符合標準	(1) IEC/EN 60601-1, IEC/EN 60601-1-2 (2) ISO 8185
Accessories / 配件	VH-3144 Rev. Humidifier Chamber w/o Scroll
Package / 包裝尺寸	12 sets per carton Carton : L65 x W46 x H36 cm