



## AUTOMATIC CONTROL HUMIDIFIER / 自控型微電腦加熱器

REF / 型號	VH-2600A-115 / VH-2600A-230
Dimensions / 規格	W135 x D170 x H156 mm
Weight / 重量	2.9 kg (w/ Chamber fitted & filled with water)
Electrical Power / 電源	(1) Supply voltage: 230±23V 115±12V (2) Supply frequency: 50/60Hz (3) Heater wire capacity: 60W max (4) Heater plate capacity: 150W, 110°C max
Heater Plate Overheat / 過熱保護	Cutoff: 115±3°C
Temperature Control / 溫度設定範圍	(1) Invasive Mode: Airway 40°C, Chamber outlet 37°C (2) Non-invasive Mode: Airway 34°C, Chamber outlet 31°C
Alarm Parameters / 警報參數	(1) Airway: High alarm / Low alarm (2) Chamber: High alarm / Low alarm
Standards / 符合標準	(1) IEC/EN 60601-1, IEC/EN 60601-1-2 (2) ISO 8185
Accessories / 配件	(1) HWP-73015 Temperature Probe (2) HWA-730F01 Heater Wire Adaptor
Package / 包裝尺寸	1 sets/ Carton L202 x W219 x H214.8 mm